

Document Code:	DPOTMH-E-57-P01-S01	
Effective Date:	06-30-2022	
Document Type:	Standard Operating Procedure	
Page Number:	1 of 4	
Department/Section:	Hematology	
Document Title:	COMPLETE BLOOD COUNT (CBC)	

PURPOSE:

To describe in detail how to process the Complete Blood Count (CBC) test correctly and always in the same manner. It is a screening procedure that is helpful in the diagnosis of many diseases.

SCOPE:

Applies to all Hematology Section Staff of Laboratory Department of Dr. Pablo O. Torre Memorial Hospital (DPOTMH)

PERSON RESPONSIBLE:

Medical Technologists (Medical Laboratory Scientists), Pathologists, Medical Doctors, Nurses, Medical Trainees, Laboratory Clerks

PROCEDURE:

- 1. Hand washing shall be performed before and after the procedure.
- 2. The wearing of Personal Protective Equipment such as gloves, masks, goggles/face shield and gowns shall be observed in doing the procedure.
- 3. The Medical Technologist collects blood in an EDTA lavender top tube (2mL/0.5mL tube) properly labeled with complete patient's name and date of collection.
- 4. The blood collected is inverted tip to tip 8-10 times gently for adequate mixing and to avoid hemolysis and destruction of chemical elements.
- Specimens collected shall be considered as biohazardous material and the disposal of syringes and needles will be according to Infection Control guidelines.
- 6. The Medical Technologist forwards properly the blood sample and patient's request to the Hematology Section. Medical Technologist then checks blood sample label, volume and clots prior to analysis. Not properly labeled sample is given back to the extractor for correct labeling. Insufficient volume and clotted



Department/Section:	Hematology	
Page Number:	2 of 4	
Document Type:	Standard Operating Procedure	
Effective Date:	06-30-2022	
Document Code:	DPOTMH-E-57-P01-S01	

sample is subject for repeat blood collection and the medical technologist, nurse incharge or the patient will be informed immediately.

7. In the HCLAB (LIS) System, the Medical Technologist enters user ID and password to display the main menu. Click Specimen Check-in icon and enter Patient Identification (PID) number (located at the upper, bottom or back part of the request) or click Find Patient Order icon and enter patient's complete name to search for the PID.

Enter PID - Date - List - Check - Accept - Print Label or Barcode

- 8. For patients with pink form or with temporary request, names are logged into a list. An ID number is assigned and written on the patient's request.
- 9. The Medical Technologist attaches the barcode sticker properly to its corresponding tube.
- 10. The Medical Technologist properly loads the blood samples on the rack & presses Start.

*for samples with 1-2mL volume – load using rack mode

*for samples with less than 1mL volume or with temporary request – use stat/manual mode

- 11. The Medical Technologist checks the abnormal results or with flaggings displayed by viewing the smear, delta checking, rerun the blood sample or repeat sample collection.
- 12. Editing is done if needed before validation. Results are auto-transferred in the HCLAB LIS then released in the BIZBOX HIS.
- 13. The Medical Technologist/Laboratory Clerk prints out and releases results to outpatients at the reception area upon presentation of the Official Receipt (OR) and results for admitted patients are can be viewed and printed out by nurses to its respective stations.



Document Code:	DPOTMH-E-57-P01-S01	
Effective Date:	06-30-2022	
Document Type:	Standard Operating Procedure	
Page Number:	3 of 4	
Department/Section:	Hematology	
Document Title:	COMPLETE BLOOD COUNT (CBC)	

- 14. Names and results of patient shall be recorded in the Hematology logbook.
- 15. Turn-around time for releasing CBC result is within 2 hours.

REFERENCE:

Hematology Handbook for Medical Technologists/Shirley I. Fabian-Cruzada, MSPH, RMT/Philippines




Document Code:	DPOTMHE-57-P01-S01	
Effective Date:	06-30-2022	
Document Type:	Standard Operating Procedure	
Page Number:	4 of 4	
Department/Section:	Hematology	
Document Title:	COMPLETE BLOOD COUNT (CBC)	

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE DIONEO-SEGOVIA, RMT Section Head, Hematology	Lads.	07-06-2022
	TIFFANY B. VILLANUEVA-COO, RMT OIC - Laboratory Manager	Elm	7-6-2022
Verified:	MONICA B. VILLANUEVA, RMT, PhD Laboratory Manager	A illanuera	7-6-2022
	MELANIE ROSE B. ZERRUDO, MD, FPSP Chair, Department of Pathology	blumafe - Jundo	7-6-2022
Reviewed:	DENNIS C. ESCALONA, MN, FPSQua Quality Assurance Supervisor	2	07-06-2022
Recommending	ROSARIO D. ABARING, MAN, PhD Ancillary Division Officer	Plaing & 1.	07.06.2022
Approval:	FREDERIC IVAN L. TING, MD OIC - Total Quality Division	B	7/8/2~
Approved:	GENESIS GOLDI D. GOLINGAN President and CEO		9/11/22



Document Title:	COMPLETE BLOOD COUNT (CBC)	
Department/Section: Hematology		
Page Number:	1 of 2	
Document Type:	Flowchart	
Effective Date:	06-30-2022	
Document Code:	DPOTMH-E-57-P01-FC01	





Document Code:	DPOTMH-E-57-P01-FC01	
Effective Date:	fective Date: 06-30-2022	
Document Type:	Flowchart	
Page Number:	2 of 2	
Department/Section:	epartment/Section: Hematology	
Document Title:	COMPLETE BLOOD COUNT (CBC)	

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE DIONEO-SEGOVIA, RMT Section Head, Hematology	Lads.	07-06-2022
	TIFFANY B. VILLANUEVA-COO, RMT OIC - Laboratory Manager	\$Un	7-6-2022
Verified:	MONICA B. VILLANUEVA, RMT, PhD Laboratory Manager	At: Uznuera	7-6-2022
	MELANIE ROSE B. ZERRUDO, MD, FPSP Chair, Department of Pathology	Dange- Generals	7-6-2022
Reviewed:	DENNIS C. ESCALONA, MN, FPSQua Quality Assurance Supervisor	0	07-06-2022
Recommending Approval:	ROSARIO D. ABARING, MAN, PhD Ancillary Division Officer FREDERIC IVAN L. TING, MD OIC - Total Quality Division	helaing k. l.	7/5/2
Approved:	GENESIS GOLDI D. GOLINGAN President and CEO		9/11/22



Document Code:	DPOTMH-E-57-P01-WI01	
Effective Date:	06-30-2022	
Document Type:	Work Instruction	
Page Number:	1 of 2	
Department/Section:	Hematology	
Document Title:	COMPLETE BLOOD COUNT (CBC)	

KEY TASKS	PERSON RESPONSIBLE
1. Encode the request of the patient system.	to BIZBOX Laboratory Clerks and Nurses
2. Extract blood sample from the patien	ts.
3. Ensure sample is acceptable for testin	g.
4. Process and analyze samples.	Medical Technologists
Release and validate results the System.	
6. Record results to Hematology logboo	ks.
7. Endorse the patient's results to the doctor.	e patient or Laboratory Clerks and Nurses



Department/Section:	Hematology	
Document Type: Page Number:	Work Instruction 2 of 2	
Effective Date: 06-30-2022		
Document Code:	nt Code: DPOTMH-E-57-P01-WI01	

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE DIONEO-SEGOVIA, RMT Section Head, Hematology	Lads.	07-06-2022
	TIFFANY B. VILLANUEVA-COO, RMT Laboratory Manager	BULW	7-6-2022
Verified:	MONICA B. VILLANUEVA, RMT, PhD Laboratory Manager	Afillznuera	7-6-2022
	MELANIE ROSE B. ZERRUDO, MD, FPSP Chair, Department of Pathology	Dange Jendo	7-6-2022
Reviewed:	DENNIS C. ESCALONA, MN, FPSQua Quality Assurance Supervisor	Ω	07-12-2022
Recommending	ROSARIO D. ABARING, MAN, PhD Ancillary Division Officer	Colaring . k. l.	07.06.2022
Approval:	FREDERIC IVAN L. TING, MD OIC - Total Quality Division	B	7/8/22
Approved:	GENESIS GOLDI D. GOLINGAN President and CEO		9/11/22