

Document Title:	PROTHROMBIN TIME (PT)
Department/Section:	Hematology
Page Number:	1 of 4
Document Type:	Standard Operating Procedure
Effective Date:	06-30-2022
Document Code:	DPOTMH-E-57-P01-S04

PURPOSE:

To describe in detail how to process the Prothrombin Time (PT) test correctly and always in the same manner. It is used to evaluate the adequacy of the extrinsic system and common pathway in the clotting mechanism and frequently used to follow the course of oral anticoagulant therapy.

SCOPE:

Applies to all Hematology Section Staff of Laboratory Department of Dr. Pablo O. Torre Memorial Hospital (DPOTMH)

PERSON RESPONSIBLE:

Medical Technologists (Medical Laboratory Scientists), Pathologists, Medical Doctors, Nurses, Medical Trainees, Laboratory Clerks

PROCEDURE:

- 1. Hand washing shall be performed before and after the procedure.
- 2. The wearing of Personal Protective Equipment such as gloves, masks, goggles/face shield and gowns shall be observed in doing the procedure.
- 3. The Medical Technologist collects blood in a 3.2% sodium citrate blue top tube (1.8mL tube) properly labeled with complete patient's name and date of collection.
- 4. The Medical Technologist collects blood sample, filled to specified volume, no underfilling or overfilling and is inverted tip to tip 3-4 times gently for adequate mixing and to avoid hemolysis and destruction of chemical elements.
- Specimens collected shall be considered as biohazardous material and the disposal of syringes and needles will be according to Infection Control guidelines.
- 6. Tests shall be performed within two hours after blood collection only.



Document Code:	DPOTMH-E-57-P01-S04
Effective Date:	06-30-2022
Document Type:	Standard Operating Procedure
Page Number:	2 of 4
Department/Section:	Hematology
Document Title:	PROTHROMBIN TIME (PT)

- 7. The Medical Technologist forwards properly the blood sample and patient's request with note of previous medications to the Hematology Section. Medical Technologist then checks blood sample label, volume and clots prior to analysis. Not properly labeled sample is given back to the extractor for correct labeling. Insufficient volume and hemolyzed/clotted sample is subject for repeat blood collection and the Medical Technologist, Nurse in-charge or the patient will be informed immediately.
- 8. Requests, unless indicated as "STAT" shall be included in the batch of examinations processed in a first come first serve basis.
- 9. Blood sample is centrifuged at a 4000 rpm for 20 minutes to ensure that the plasma is platelet poor.
- 10. In the HCLAB LIS system, the Medical Technologist enters user ID and password to display the main menu. Click Specimen Check-in icon and enter Patient Identification (PID) number (located at the bottom, upper or back part of the request) or click Find Patient Order icon and enter patient's complete name to search for the PID.

Enter PID - Date - List - Check - Accept - Print Label or Barcode

- 11. The Medical Technologist logs into a list the patient names. An ID number is assigned and written on the patient's request.
- 12. The Medical Technologist attaches the barcode sticker properly to its corresponding tube.
- 13. The Medical Technologist properly loads the blood samples on the rack & presses Start.
- 14. Results of all parameters are displayed and copied in the machine's monitor.



Document Code:	DPOTMH-E-57-P01-S04	
Effective Date:	06-30-2022	
Document Type:	Standard Operating Procedure	
Page Number:	3 of 4	
Department/Section:	Hematology	
Document Title:	PROTHROMBIN TIME (PT)	

- 15. The Medical Technologist checks the abnormal results or with flaggings displayed by viewing the smear, delta checking or rerun the blood sample or correlate result to patient's anticoagulant therapy.
- 16. Editing is done if needed before validation. Results are auto transferred in the HCLAB LIS then released in the BIZBOX HIS system.
- 17. The Medical Technologist/Laboratory Clerk prints out and releases results to outpatients at the reception area upon presentation of the Official Receipt (OR) and results for admitted patients are can be viewed and printed out by nurses to its respective stations.
- 18. Critical PT results must be reported to the Attending Physician or Nurse immediately.
- 19. Names and results of patient shall be recorded in the Hematology logbook.
- 20. Turn-around time for releasing PT result is within 2 hours.

REFERENCE:

Hematology Handbook for Medical Technologists/Shirley I. Fabian-Cruzada, MSPH, RMT/Philippines



Document Code:	DPOTMH-E-57-P01-S04
Effective Date:	06-30-2022
Document Type:	Standard Operating Procedure
Page Number:	4 of 4
Department/Section:	Hematology
Document Title:	PROTHROMBIN TIME (PT)

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE DIONEO-SEGOVIA, RMT Section Head, Hematology	clads.	07-06-2022
	TIFFANY B. VILLANUEVA-COO, RMT Laboratory Manager	Pth	7-6-2022
Verified:	MONICA B. VILLANUEVA, RMT, PhD Laboratory Manager	D'llarun	7-6-2022
	MELANIE ROSE B. ZERRUDO, MD, FPSP Chair, Department of Pathology	bluenge-Gensk	7-4-2022
Reviewed:	DENNIS C. ESCALONA, MN, FPSQua Quality Assurance Supervisor	2	07-14-2022
Recommending Approval:	ROSARIO D. ABARING, MAN, PhD Ancillary Division Officer FREDERIC IVAN L. TING, MD	flaing b. 4.	07.06.2022
Approved:	OIC - Total Quality Division GENESIS GOLDI D. GOLINGAN President and CEO		9/11/22



Document Code:	DPOTMH-E-57-P01-WI04
Effective Date:	06-30-2022
Document Type:	Work Instruction
Page Number:	1 of 2
Department/Section:	Hematology
Document Title:	PROTHROMBIN TIME (PT)

KEY TASKS		PERSON RESPONSIBLE
1.	Encodes the request of the patient to BIZBOX system.	Laboratory Clerks and Nurses
2.	Extracts blood sample from the patients.	
3.	Ensures sample is acceptable for testing.	
4.	Disposes properly and safely biohazardous and infectious wastes and materials.	
5.	Ensures laboratory equipment is clean, functional and calibrated.	Medical Technologist/Nurse
6.	Processes and analyzes samples. Reports result with significant findings.	
7.	Releases and validates result thru BIZBOX System.	
8.	Records result to Hematology logbooks.	
9.	Endorses the patient's results to the patient or doctor.	Laboratory Clerks and Nurses




Document Title:	PROTHROMBIN TIME (PT)
Department/Section:	Hematology
Page Number:	2 of 2
Document Type:	Work Instruction
Effective Date:	06-30-2022
Document Code:	DPOTMH-E-57-P01-WI04

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE DIONEO-SEGOVIA, RMT Section Head, Hematology	afads.	07-06-2022
	TIFFANY B. VILLANUEVA-COO, RMT Laboratory Manager	toller	7-6-2022
Verified:	MONICA B. VILLANUEVA, RMT, PhD Laboratory Manager	D: blancer	7-6-2022
	MELANIE ROSE B. ZERRUDO, MD, FPSP Chair, Department of Pathology	blenge- Gend	7-6-2022
Reviewed:	DENNIS C. ESCALONA, MN, FPSQua Quality Assurance Supervisor	2	07-04-2022
Recommending Approval:	ROSARIO D. ABARING, MAN, PhD Ancillary Division Officer FREDERIC IVAN L. TING, MD	Plainz & 1.	07.06.2022
Approved:	OIC - Total Quality Division GENESIS GOLDI D. GOLINGAN President and CEO		9/11/22



Document Code:	DPOTMH-E-57-P01-FC04	
Effective Date:	06-30-2022	
Document Type:	Flowchart	
Page Number:	1 of 2	
Department/Section:	Hematology	
Document Title:	PROTHROMBIN TIME (PT)	





Document Code:	DPOTMH-E-57-P01-FC04
Effective Date:	06-30-2022
Document Type:	Flowchart
Page Number:	2 of 2
Department/Section:	Hematology
Document Title:	PROTHROMBIN TIME (PT)

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE DIONEO-SEGOVIA, RMT Section Head, Hematology	ofads.	07-06-2022
	TIFFANY B. VILLANUEVA-COO, RMT Laboratory Manager	Blw	7-6-2022
Verified:	MONICA B. VILLANUEVA, RMT, PhD Laboratory Manager	Aillanuer-	7-6-2022
	MELANIE ROSE B. ZERRUDO, MD, FPSP Chair, Department of Pathology	bluerge - Generale	7-6-2022
Reviewed:	DENNIS C. ESCALONA, MN, FPSQua Quality Assurance Supervisor	0	07-06-0022
Recommending Approval:	ROSARIO D. ABARING, MAN, PhD Ancillary Division Officer FREDERIC IVAN L. TING, MD	Maing & l.	07.06.2022
Approvai.	OIC - Total Quality Division	1	7/8/22
Approved:	GENESIS GOLDI D. GOLINGAN President and CEO	-	9/11/22