43		THE REAL PROPERTY.
3.		7
	1	

Document Code:	DPOTMH-E-60-P05-S01
Effective Date:	08-30-2022
Document Type:	Standard Operating Procedure
Page Number:	1 of 3
Department/Section:	Molecular Laboratory
Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS

PURPOSE:

To standardize practices in washing benches and floors while maintaining a clean environment, thus, reducing infectious risk to Laboratory Staff as well as reducing the risk of cross-contamination.

SCOPE:

All Molecular Laboratory Personnel of Dr. Pablo O. Torre Memorial Hospital (DPOTMH)

PERSON RESPONSIBLE:

Medical Technologist, Laboratory Aide, Housekeeping Personnel

PROCEDURE:

- 1 Always use the appropriate protective gear.
- 2 Laboratory workbenches must be cleaned between uses. Keeping the laboratory bench clean ensures leftover dirt or other materials do not contaminate future processes. It can also ensure that dangerous substances do not stay longer on surfaces that may cause damage to the equipment or cause harm to other personnel.
- 3 Laboratory floors should be mopped at least daily. In doing so, first make use of this:
 - 3.1 Aisles are clear of boxes, supplies, or other obstructions.
 - 3.2 Loose wires, cables, and computer cords are tied and organized.
- 4 Emergency areas like eyewash stations, showers, and fire extinguishers should remain unobstructed at all times.
- 5 Laboratory areas should be dusted and decluttered regularly.
- 6 Clean items that may not be part of regular testing or experiments. Wipe down chairs, telephones, computers, timers, pens, etc. daily to make sure contagions aren't accidentally left behind.



Document Code:	DPOTMH-E-60-P05-S01
Effective Date:	08-30-2022
Document Type:	Standard Operating Procedure
Page Number:	2 of 3
Department/Section:	Molecular Laboratory
Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS

- 7 Remove loose items from the Laboratory Work Station or cabinets that you are working with to eliminate the possibility of missing concealed areas during disinfection. Cold racks, pipettes, etc. should be relocated and washed appropriately.
- 8 All corners, edges, and undersides should be cleaned.
- 9 Label diluted cleaning solutions. Store properly.
- 10 Use 70% ethanol while cleaning equipment surfaces. Floors can be mopped using a hypochlorite solution.
- All personnel must practice CLAYGO (Clean as you Go) and must observe proper waste decontamination and disposal.

REFERENCES:

- List N: Disinfectants for use against SARs-CoV-2, https://cfpub.epa.gov/giwiz/disinfectants/index.cfm
- 2. https://www.nationalhogfarmer.com/health-diseases/1015-effective-cleaning-disinfecting-steps
- 3. https://www.cdc.gov/coronavirus/2019-ncov/community/disinfecting-building-facility.html
- 4. https://www.labtechsupplyco.com/how-to-properly-clean-and-disinfect-your-lab-benches/



Document Code:	DPOTMH-E-60-P05-S01	
Effective Date:	08-30-2022	
Document Type:	Standard Operating Procedure	
Page Number:	3 of 3	
Department/Section:	Molecular Laboratory	
Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS	

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE D. SEGOVIA, RMT	afads.	09-05-2022
	Section Head, Hematology	C	
	MONICA B. VILLANUEVA, RN, RMT, PhD	IC.II	0.5.00
	Laboratory Director	Willemen	77-2022
Verified by:	DAVID G. PEDROZA, MD	Samuel	95-22
	Associate Pathologist		-
	MELANIE ROSE B. ZERRUDO, MD, FPSP	lever for fout	9-5-102
	Chair, Department of Pathology	18	1
Reviewed by:	DENNIS C. ESCALONA, MN, FPCHA, FPSQua	0	9-05-2022
	Quality Assurance Supervisor		
	ROSARIO D. ABARING, RN, MN, PhD, FPCHA	01.1	1
Recommending	Ancillary Services Division Officer	Il Warmer &	.09.05.20
Approval by:	MA. ANTONIA S. GENSOLI, MD, FPPS, FPCHA	"	@ 22 2Y
	Vice President – Chief Medical Officer	many	9.27.27
Approved by:	GENESIS GOLDI D. GOLINGAN	HAD.	Palaelas
11	President and CEO	100	11200



Document Code:	DPOTMH-E-60-P05-WI01
Effective Date:	08-30-2022
Document Type:	Work Instruction
Page Number:	1 of 2
Department/Section:	Molecular Laboratory
Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS

KEY TASKS	PERSON RESPONSIBLE
1. Cleans the laboratory workbenches between uses.	
2. Ensures that dangerous substances do not stay longer on surfaces that may cause damage to the equipment or cause harm to other personnel.	
3. Mops laboratory floors at least daily.	
4. Clears aisles of boxes, supplies, or other obstructions.	Housekeeping Personnel
5. Ensures emergency areas like eyewash stations, showers, and fire extinguishers are unobstructed.	
6. Cleans items that may not be part of regular testing or experiments.	
7. Removes loose items from the Laboratory Work Station or cabinets to omit the possibility of missing concealed areas during disinfection.	
8. Cleans all corners, edges, and undersides.	




Document Code:	DPOTMH-E-60-P05-WI01	
Effective Date:	08-30-2022	
Document Type:	Work Instruction	
Page Number:	2 of 2	
Department/Section:	Molecular Laboratory	
Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS	

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE D. SEGOVIA, RMT	-fads ·	09-05-2022
	Section Head, Hematology	(-	
	MONICA B. VILLANUEVA, RN, RMT, PhD	11.24	
	Laboratory Director	Mr. Unner.	9-5-30-31
Verified by:	DAVID G. PEDROZA, MD	Sulley	650
	Associate Pathologist	ENUI)	
	MELANIE ROSE B. ZERRUDO, MD, FPSP	Quele Cars	19-5-202
	Chair, Department of Pathology	18	,
Reviewed by:	DENNIS C. ESCALONA, MN, FPCHA, FPSQua	0	9-22-2022
	Quality Assurance Supervisor		
	ROSARIO D. ABARING, RN, MN, PhD, FPCHA	M. 11	
Recommending	Ancillary Services Division Officer	Hylaing & 4	.09.05.20
Approval by:	MA. ANTONIA S. GENSOLI, MD, FPPS, FPCHA		
11	Vice President – Chief Medical Officer	man	9-27.21
Approved by:	GENESIS GOLDI D. GOLINGAN	And Co	dal
••	President and CEO	100	7/20/22




Document Code:	DPOTMH-E-60-P05-FC01	
Effective Date:	08-30-2022	
Document Type:	Flowchart	
Page Number:	1 of 2	
Department/Section:	Molecular Laboratory	
Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS	





Document Code:	DPOTMH-E-60-P05-FC01	
Effective Date:	08-30-2022	
Document Type:	Flowchart	
Page Number:	1 of 2	
Department/Section:	Molecular Laboratory	
Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS	





Document Title:	CLEANING AND DISINFECTING OF MOLECULAR BIOLOGY LABORATORY WORK AREAS
Department/Section:	Molecular Laboratory
Page Number:	2 of 2
Document Type:	Flowchart
Effective Date:	08-30-2022
Document Code:	DPOTMH-E-60-P05-FC01

APPROVAL:

	Name/Title	Signature	Date
Prepared by:	REDABELLE D. SEGOVIA, RMT	dads.	09-05-2022
	Section Head, Hematology		
Verified:	MONICA B. VILLANUEVA, RN, RMT, PhD	drillamen	9-5-22
	Laboratory Director		
	DAVID G. PEDROZA, MD	SARTH	9/5/22
	Associate Pathologist		
	MELANIE ROSE B. ZERRUDO, MD, FPSP	Shen fe fens	lo 9-5-22
	Chair, Department of Pathology		
Reviewed by:	DENNIS C. ESCALONA, MN, FPCHA, FPSQua	0	9-22-2002
	Quality Assurance Supervisor		
	ROSARIO D. ABARING, RN, MN, PhD, FPCHA	1.	,
Recommending	Ancillary Services Division Officer	Helains b.	. 07.05.20
Approval:	MA. ANTONIA S. GENSOLI, MD, FPPS, FPCHA		
	Vice President- Chief Medical Officer	make	9-27-
Approved:	GENESIS GOLDI D. GOLINGAN	And the	Palaelan
	President and CEO	TONIO	1/2/02